

---

Poster

## [P27-5] P27-5: Cardiovascular drugs (2)

Chair: David A. Joyce, Australia

Wed. Sep 27, 2017 12:30 PM - 1:30 PM Annex Hall (1F)

---

(Wed. Sep 27, 2017 12:30 PM - 1:30 PM Annex Hall )

## [P27-5-7] Determination of homocysteine in human plasma / serum using LC-MS8060

Changkun Li (Shimadzu (China) CO., LTD.)

Keywords: Homocysteine, Plasma, LCMS-8060

### Background

Homocysteine (Hcy) is the intermediate for methionine metabolism. It is biosynthesized from methionine and a homologue of the amino acid cysteine, differing by an additional methylene bridge (-CH<sub>2</sub>-). An enhanced blood homocysteine concentration is a medical condition which is termed as hyperhomocysteinemia. It is considered as an independent risk factor for degenerative vascular diseases, which may lead to heart attack, stroke and thromboses.

### Methods

A combined system of Shimadzu LC-30A and triple quadrupole mass spectrometer LCMS-8060 was used in the experiment. Homocysteine was measured using the commercially available ClinMass® Complete Kit for total homocysteine in Plasma and Serum. Calibrators, control samples were provided by the kit. The samples are spiked with an internal standard (50  $\mu$ l), the bound homocysteine is released (reduction) and after a precipitation step. The mass spectrometer was operated in the positive mode. Quantification was performed using MRM-Mode of the transitions of m/z 136.2 - m/z 90.2. The mobile phase was acetonitrile(0.1% formic) and water(0.1% formic), a Discovery HS F5-3 (2.1 mmI.D.  $\times$ 150 mmL, 3 m) column was used for analysis samples.

### Results

Using the LCMS-8060 analysis of Homocysteine in Plasma and Serum can be performed quite easily. The calibration curves showed good linearity within the range of 0.794~6.86 mg/L with the correlation coefficient not less than 0.9996. The Method Detection Limit (MDL) and Method Quantification Limit (MQL) is 0.26 ng/mL and 0.87 ng/mL respectively, and the accuracy of quality control samples are in the range of target value.

### Conclusions

A selective and sensitive ultra-high performance liquid chromatography - triple quadrupole mass spectrometry (UHPLC-MS/MS) method was developed for the direct determination of homocysteine in human plasma/serum. This method showed good sensitivity and linearity in a clinically relevant concentration range, and suitable for detection of clinical samples.