
Poster

[P27-8] P27-8: Assay and monitoring

Chair: Yoshihiko Hirotsu, Japan

Wed. Sep 27, 2017 12:30 PM - 1:30 PM Annex Hall (1F)

(Wed. Sep 27, 2017 12:30 PM - 1:30 PM Annex Hall)

[P27-8-7] A new LC-MS/MS method for determination of vitamin K1 and vitamin K2 (MK-4, MK-7) in serum

Eva Klappkova¹, Jana Cepova², Katerina Dunovska³, Milada Halacova⁴, Richard Prusa⁵ (1.Charles University and University Hospital Motol, 2.Charles University and University Hospital Motol, 3.Charles University and University Hospital Motol, 4.Charles University and University Hospital Motol, 5.Charles University and University Hospital Motol)

Keywords: vitamin K1, vitamin K2, LC-MS

Background

Vitamins K1 and K2 are the cofactors for the enzyme gamma-glutamylcarboxylase, which is involved in carboxylation of the vitamin K-dependent proteins. The role of vitamin K extends well beyond the regulation of blood clotting to impact bone formation, and development of heart disease, and possibly cancer.

Objective of this study was to evaluate a new LC-MS/MS method for determination of vitamin K1 and two forms of vitamin K2 (MK-4, MK-7) in patient serum.

Methods

A UPLC-tandem mass spectrometry method for the determination of vitamin K1 and vitamin K2 (MK-4 and MK-7) in human serum was developed. Mass spectrometric detection was performed on a triple-quadrupole in the positive electrospray ionization mode by multiple reaction monitoring (Agilent 1290 with Triple Quad 6470, Agilent Technologies). K1-d7 and MK-4-d7 were used as internal standards. 20 μ L of IS mixture were added to 500 μ L of serum and 2 mL of ethanol were added to precipitate the proteins. The mixture was extracted with 4 mL of hexane for 10 min and then centrifuged. Solid phase extraction was then used. Chromatographic separation was carried out on SB-C8 column (2.1x100 mm, 1.8 μ m, Agilent Technologies) using a binary gradient of mobile phases (A - 0.1% ammonium fluoride in 50% methanol, B - 0.1% ammonium fluoride in methanol). The separation was accomplished at 40 °C. The analysis time was 9 minutes.

Results

The LC-MS/MS method has been successfully validated. The linearity was excellent up to 10 ng/mL for all three forms of vitamin K. Lower limits of quantitation were 0.03 ng/mL for K1 and MK-4, and 0.05 ng/mL for MK-7. The intra and interday accuracy and precision were evaluated on two QC samples by multiple analysis and coefficients of variation were 2.3-7.0% for intra-assay, 7.2-11.3% for inter-assay, and 0.5-6.6% for precision.

Conclusions

We developed and fully validated a new LC-MS/MS method for determination of vitamin K1, MK-4 and MK-7, which could be implemented into the practice, because new studies have reported that daily vitamin K intake improves bone metabolism, vascular calcification, and cardiovascular risk, and this method allows the determination of the optimal dose.

Supported by the project (Ministry of Health, Czech Republic) for conceptual development of research organization 00064203 (University Hospital Motol, Prague, Czech Republic)

